



# Opticyte<sup>®</sup>

Amniotic Ocular **Matrix**


More than just a membrane, Opticyte Ocular Matrix is designed to preserve the native amnion with a basement membrane similar to the conjunctiva.<sup>1</sup>

Opticyte Ocular Matrix is derived from human amnion that has been gently processed and sterilized using an ISO 11137 validated terminal sterilization with a low dose of irradiation.

This customized process yields a graft that is optimized for use by ophthalmic professionals as a biological bandage that protects and provides a protein-rich extracellular matrix. Indication is for homologous use for covering and protecting:

- Corneal and stromal ulcerations
- Pterygium excisions
- Conjunctiva surface reconstruction<sup>2</sup>

For this reason, current theory suggests that amniotic membrane augments support for epithelial cells, limbal stem cells, and corneal transient amplifying cells.<sup>3</sup>



Processed Amnion Histology and Microscopy.

Source data on file

## PROCEDURE

Opticyte Ocular Matrix can be easily applied to the ocular surface by:

- Absorbable sutures
- Biologic tissue adhesive
- Contact lenses<sup>4</sup>

## SAFETY & QUALITY

Opticyte Ocular Matrix is processed at an FDA-registered and AATB-accredited institution, providing clinicians a safe and high-quality tissue product for use in ophthalmic procedures.

For more information related to Opticyte insurance benefit verifications, prior authorization assistance, or claims appeal assistance, please contact our **Reimbursement Support Line**.

**Phone:** 919.921.8105

**Fax:** 919.267.3753

**Email:** support@merakris.com

## Clinical

## Surgical

Item Number:	Size	Item Number:	Size
AM-OP1008	8 mm	AM-OPS1011	1x1 cm <sup>2</sup>
AM-OP1010	10 mm	AM-OPS1012	1x2 cm <sup>2</sup>
AM-OP1012	12 mm	AM-OPS1014	1x4 cm <sup>2</sup>
AM-OP1014	14 mm	AM-OPS1023	2x3 cm <sup>2</sup>
AM-OP1016	16 mm		

### Opticyte Graft Placement Instruments

OPI-FCPS-01	Tissue Forceps Stainless Steel
OPI-SPEC-01	Speculum
OPI-STRY	Instrument Sterilization Caddy



## Quality Compliance & Donor Screening

Opticyte tissue is procured and processed in the United States according to standards and/or regulations established by the American Association of Tissue Banks (AATB) and the United States Food & Drug Administration (FDA). Donor ethics and negative FDA-approved serological screening includes:

- Hepatitis B Surface Antigen (HBsAg)
- Hepatitis B core Antibody (HBcAb)
- Hepatitis C Antibody (HCV)
- HIV 1/2 Plus O Antibody
- Syphilis Serological Detection Test
- HTLV 1-2 Antibody
- HCV NAT (Nucleic Acid Test)
- HBV NAT (Nucleic Acid Test)
- WNV NAT (Nucleic Acid Test)

### References

1. Fukuda K et al. Differential distribution of subchains of the basement membrane components type IV collagen and laminin among the amniotic membrane, cornea, and conjunctiva. *Cornea*. 199;18:73-79.
2. Tseng SCG et al. Amniotic membrane transplantation with or without limbal allografts for cornea surface reconstruction in patients with limbal stem cell deficiency. *Arch Ophthalmol*. 1998;116:431-441.
3. "Amniotic Membrane Transplantation" EyeWiki, American Academy of Ophthalmology, December 19, 2017, [http://eyewiki.aao.org/Amniotic\\_Membrane\\_Transplant#Procedure](http://eyewiki.aao.org/Amniotic_Membrane_Transplant#Procedure)
4. Tseng SCG et al. *Amniotic membrane suturing techniques*. 2007. Springer Berlin Heidelberg.

### Regulatory Factors

Opticyte Ocular Matrix is classified and qualifies as a human tissue allograft (HCT/P) as outlined in 21CFR 1271 under Section 361 of the Public Health Service Act. Opticyte Ocular Matrix is intended for "Homologous Use" as it is used to cover and protect tissue.

### Disclaimer

This product is intended for homologous use and should not be used to prevent, treat, or cure a disease. Any use of this product outside of homologous use area may require an Investigational New Drug or Investigational Device Exemption filing with the FDA.